



MICROARRAY SPOTTING PIN OFFERINGS

Available for purchase worldwide (except the U.S.)

All pins licensed under Incyte US patent Nos.
 #5,807,522, #6,110,426, EP #0 804 731, AU #199528629, CA #2,192,095 and JP #10503841

All pins are manufactured in Point Technologies' new state-of-the-art facility located in Global Park, Heredia, Costa Rica.



ACCELERATOR™ brand Tungsten microarray **split** pins in PT and PTL configurations:

PT overall length 1.061" (26.95 mm)

PTL overall length 1.75" (44.45 mm)

<u>CATALOG NUMBER</u>	<u>PIN DIAMETER</u>	<u>TIP DIAMETER</u>	<u>SPOT DIAMETER</u>	<u>END USER PRICE</u>
PT/PTL 2000	0.061" (1.55 mm)	65µm (0.0025")	67-90um	\$225.00
PT/PTL 2500	0.061" (1.55 mm)	74µm (0.0029")	75-100um	\$225.00
PT/PTL 3000	0.061" (1.55 mm)	89µm (0.0035")	90-120um	\$225.00
PT/PTL 4000	0.061" (1.55 mm)	119µm (0.0046")	120-155um	\$225.00
PT/PTL 5000	0.061" (1.55 mm)	147µm (0.0057")	150-215um	\$225.00
PT/PTL 6000	0.061" (1.55 mm)	175µm (0.0069")	180-215um	\$225.00



ACCELERATOR™ brand Tungsten microarray **solid** pins in PTS and PTLS configurations:

PTS overall length 1.061" (26.95 mm)

PTLS overall length 1.75" (44.45 mm)

<u>CATALOG NUMBER</u>	<u>PIN DIAMETER</u>	<u>TIP DIAMETER</u>	<u>SPOT DIAMETER</u>	<u>END USER PRICE</u>
PTS/PTLS 50	0.061" (1.55 mm)	50µm (0.0020")	75-100um	\$99.00
PTS/PTLS 100	0.061" (1.55 mm)	100µm (0.0039")	100-125um	\$99.00
PTS/PTLS 200	0.061" (1.55 mm)	200µm (0.0079")	140-150um	\$99.00
PTS/PTLS 300	0.061" (1.55 mm)	300µm (0.0118")	220-250um	\$99.00
PTS/PTLS 400	0.061" (1.55 mm)	400µm (0.0158")	280-300um	\$99.00
PTS/PTLS 500	0.061" (1.55 mm)	500µm (0.0197")	360-400um	\$99.00
PTS/PTLS 600	0.061" (1.55 mm)	600µm (0.0236")	430-450um	\$99.00

We look forward to working with you.

P 800.557.7059
 sales@pointtech.com

f 303.415.9866
 www.pointtech.com

Point Technologies, Inc.
 6859 N. Foothills Hwy.
 Boulder, CO 80302
 USA



PFT™ Stainless Steel Microarray Pins for Stealth-style pin holders.

Dimensions: Length = 1.061" (26.95 mm), Dia = 0.061" (1.55 mm).

Also available in Tungsten.

<u>CATALOG NUMBER</u>	<u>TIP DIAMETER</u>	<u>SPOT DIAMETER</u>	<u>END USER PRICE</u>
PFT-2500	85µm (0.0035")	85-100 µm	\$199.00
PFT-3000	100µm (0.0039")	100-125 µm	\$199.00
PFT-4000	110µm (0.0043")	135-150 µm	\$199.00



PCH™ Stainless Steel Microarray Pins for Chipmaker™-style pin holders.

Dimensions: Length = 1.750" (44.45 mm), Dia = 0.061" (1.55 mm)

<u>CATALOG NUMBER</u>	<u>TIP DIAMETER</u>	<u>SPOT DIAMETER</u>	<u>END USER PRICE</u>
PCH-3000	100µm (0.0039")	100-125 µm	\$199.00
PCH-4000	110µm (0.0043")	135-150 µm	\$199.00



Point Technologies PMJ™ slit pins for Brown Lab-style pin holders made from 440C SS.

Dimensions: Length = 2.50" (63.5 mm), Dia = 0.061" (1.55 mm)

<u>CATALOG NUMBER</u>	<u>TIP DIAMETER</u>	<u>SPOT DIAMETER</u>	<u>END USER PRICE</u>
PMJ-2000	75 -600 µm per request	75 -800 µm per request	\$199.00

Chipmaker™ - Trade marked by TeleChem International, Inc.